




DEGASI® PLUS CLASSIC compared to the old DEGASI® CLASSIC

DEGASI® PLUS SERIES

Our latest addition to the DEGASI® standalone series is the DEGASI® PLUS. It can be equipped for a wide range of applications, depending on what kind of performance you need.

"It's small, smart & silent..."



THE FIRST CHOICE

DEGASI® PLUS is a further development of degassers such as:

- DEGASI® Classic
- DEGASI® Micro
- DEGASI® GPC
- DEGASI® Semi-Prep

FEATURES

- Much smaller footprint, height and weight
- New vacuum pump using a smart board
- Allmost silent
- Available with 1-6 channels
- Parallel connected vacuum chambers
- Easier accessibility and visibility
- Easier troubleshooting and replacement of parts
- Short leadtime

The old DEGASI® CLASSIC and old DEGASI® GPC are still available:



Old DEGASI® CLASSIC

Old DEGASI® GPC

Part Number	Number of Channels	Internal Volume	Part number	Number of Channels	Internal Volume
0001-6352-A	2	480 µl	0001-6622	2	480 µl
0001-6353-A	3	480 µl	0001-6623	3	480 µl
0001-6354-A	4	480 µl	0001-6624	4	480 µl
0001-6355-A	5	480 µl	0001-6625	5	480 µl



All DEGASI® PLUS Units are delivered with a Systec ZHCR Vacuum Pump with Smart Board